

New Approach Methodologies (NAMs)

1 in 3
MAbs tested
are polyspecific




We identified off-target binding in 33% of lead candidate MAbs and in 15% of clinically approved antibodies, highlighting a critical need for improved specificity testing

Scientifically validated NAMs provide an opportunity to better predict drug safety in humans while also accelerating drug development and regulatory review

Norden et al., 2024, mAbs

Membrane Proteome Array™ (MPA)


Better Specificity Screening Can De-Risk Antibody Development



- Compatible with all biotherapeutic modalities
- Proteins screened under unfixed conditions to maintain native conformations
- MPA screening enables confident lead selection and better specificity data for IND
- Optional additional 1,200 native secreted proteins screened by ELISA


Assay Validation for DDT Qualification

1. Library Based on FDA-Recommended Tissues for Specificity Testing



- All targets are tested for binding while expressed within human cells
- Proteins screened under unfixed conditions to maintain native conformations

2. Quantitative Data to Analyze Results and Justify Interpretations



3. Analytical Validation for Full Qualification Package (FQP)

Reproducibility


Intermediate Precision Analysis
• Across days, concentrations, lots, operators

Internal Controls Over Time
• 2 MAbs tested repeatedly over years

Quality control


Assay Validation
Specific target identification:
• False-negative rate <1% calculated in 2 different ways
• False-positive hits eliminated in validation

Sensitivity:
• Conditions optimized for high sensitivity



Clinical Validation: Comparison to Tissue Cross-Reactivity

1. MPA Identified Off-Target Binding in Clinical Antibodies



Norden et al., 2024, mAbs

- Clinical-stage and FDA-approved MAbs were tested on the MPA
- The MPA identified polyspecificity at all stages of clinical development, underscoring the critical need for improved specificity testing

2. MPA Superior Over TCR in Identifying Off-Targets

# MAbs	MPA result	TCR result
7	Off-target binding	Clean
3	Clean	Unexpected cross-reactivity
17	Clean	Clean
8	Clean	Inconclusive (additional cytoplasmic staining, technical limitations, variability, false-positive staining)

- An unbiased selection of 35 FDA-approved antibody therapeutics were screened on the MPA and results were compared to TCR results obtained from public BLA reviews
- Ambiguous TCR results highlight technical limitations and the lack of actionable data resulting from many TCR studies

MPA identified off-targets missed in TCR studies. Preliminary data suggests MPA results are better correlated with clinical toxicity

MPA in Final Review as an FDA-Qualified Drug Development Tool

- FDA guidance recommends the use of advanced, human-relevant, in vitro assays (NAMs) to assess antibody specificity and drug safety
- MPA technology in final stages of DDT qualification (FQP submitted)
- MPA's ongoing process enhancements incorporate direct feedback from the FDA through its ISTAND program



Integral Molecular is the industry leader in working with membrane proteins, backed by 20+ years of experience, 600+ customers, and 500+ manuscripts using our technologies.

MPA services are available directly through Integral Molecular or via our partners and international distributors.

Why Choose the MPA

- TRUSTED** – Successfully screened 2,000+ MAbs using ISO 9001 certified quality processes; trusted by hundreds of customers
- BIOLOGICALLY RELEVANT** – The only screening platform that assesses its entire library under native conditions
- ACCEPTED** – Data accepted by FDA, EMA, NMPA; 100+ IND submissions filed with MPA data
- QUANTITATIVE** – Flow cytometry provides sensitive and quantitative measurements

v1025

