


### Off-Target Binding Poses Safety Risks



**1 in 3**


MAbs tested are polyspecific

Off-target binding can cause:

- Unintended cellular toxicity
- Adverse events
- Clinical trial failures

### Membrane Proteome Array™ (MPA)

*A Better Specificity Test To De-Risk Antibody Development*



- Cell-based, native-conformation protein array
- Comprehensive human membrane proteome
- High-throughput assessment of antibody specificity using quantitative flow cytometry

### MPA in Final Review to Become an FDA-Qualified NAM

Validation process for qualification through FDA's IStand program

#### Analytical Validation

Validation of platform, protocols, procedures and justification of data analysis

- Performance characteristics (target identification rate, false-negative rate)
- Precision, sensitivity, reproducibility

#### Clinical Validation

Evaluate overall utility for clinical relevance and drug development


- Prospective and retrospective studies on pre-clinical and clinical samples
- Comparison to current standard



\*MPA's Full Qualification Package (FQP) currently under review

## A Quantitative, Human-Relevant NAM to Test Specificity

### 1. Cell-Based Protein Array of ~6,000 Native Human Membrane Proteins



Membrane Proteome Array:  
5,220 unique proteins


94% coverage  
of membrane proteome

#### Key features:

- ✓ V5 tag to confirm expression
- ✓ Heterocomplexes
- ✓ GPI-linked proteins
- ✓ Sex-specific proteins included
- ✓ Developmental stages represented
  - Fetal proteins
  - Placental proteins

- All targets are tested for binding while expressed within human cells
- Additional 1,200 native secreted proteins screened by ELISA

### 2. Quantitative Binding Data




- Compatible with all biotherapeutic modalities

### 3. Specificity Data Routinely Accepted by FDA for IND Submissions

- MBRC-101 is a first-in-class antibody-drug conjugate (ADC) targeting EphA5, a receptor tyrosine kinase expressed in multiple solid tumors
- Anti-EphA5 antibody tested on MPA
- FDA granted IND, no tissue cross-reactivity study required



MBRC-101 currently in Phase 1/2 (NCT06014658)




Staquinini et al., 2025, J Clin Invest

MPA provides confidence for lead selection and specificity data for IND

## Comparison to Tissue Cross-Reactivity (TCR) Studies

### 1. MPA Identified Off-Target Binding in Clinical Antibodies



#### Clinical MABs Screened on the MPA

|            | All | Approved | Phase 2/3 | Withdrawn |
|------------|-----|----------|-----------|-----------|
| # MABs     | 83  | 40       | 25        | 18        |
| Off-target | 15  | 6        | 5         | 4         |
| Rate       | 18% | 15%      | 20%       | 22%       |

Norden et al., 2024, mAbs

- Clinical-stage and FDA-approved MABs were tested on the MPA
- The MPA identified polyspecificity at all stages of clinical development, underscoring the critical need for improved specificity testing

### 2. MPA Superior Over TCR in Identifying Off-Target Safety Risks

| # MABs | MPA result                                                                                                                         | TCR result                                                                                                     |
|--------|------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| 7      | Off-target binding                                                                                                                 | Clean                                                                                                          |
|        | 4 have FDA black box warnings<br>2 have FDA warnings of severe side-effects<br>1 shows no toxicity, but MPA off-target is in brain |                                                                                                                |
| 3      | Clean                                                                                                                              | Unexpected cross-reactivity                                                                                    |
|        | In all 3 cases no off-target toxicity was identified in NHP or human studies                                                       |                                                                                                                |
| 17     | Clean                                                                                                                              | Clean                                                                                                          |
| 8      | Clean                                                                                                                              | Inconclusive<br>(additional cytoplasmic staining, technical limitations, variability, false-positive staining) |


- An unbiased selection of 35 FDA-approved antibody therapeutics were screened on the MPA and results were compared to TCR results obtained from public BLA reviews
- Ambiguous TCR results highlight technical limitations and the lack of actionable data associated with many TCR studies

MPA identified off-targets missed in TCR studies; preliminary data suggests MPA results are better correlated with clinical toxicity

## Why Choose the MPA




**TRUSTED** – Successfully screened 2,500+ MABs using ISO 9001 certified quality processes; trusted by hundreds of customers



**BIOLOGICALLY RELEVANT** – The only screening platform that assesses its entire library under native conditions



**ACCEPTED** – Data accepted by FDA, EMA, NMPA; 100+ IND submissions filed with MPA data



**QUANTITATIVE** – Flow cytometry provides sensitive and quantitative measurements

v0426

