

Infection of 293T-hsACE2 cells with SARS-CoV-2 firefly luciferase Reporter Virus Particles

Safety note: SARS-CoV-2 RVPs are created using a second-generation lentiviral system and should be used under BSL2 or enhanced BSL2 laboratory conditions. Use of SARS-CoV-2 RVPs requires infectious agent training. Always wear gloves, eye protection, and a lab coat when using RVPs. Filter tips on pipettors are recommended when handling infectious agents. Please see the SARS-CoV-2 RVP Safety Data Sheet for full safety information. <https://www.integralmolecular.com/sarscov2rvpddocuments/>

Required Reagents

- SARS-CoV-2 firefly luciferase Reporter Virus Particles (RVPs) (Integral Cat# RVP-7xxF)
- Clonal 293T-hsACE2 cells (Integral Cat# C-HA102)
- 96 or 384 well plate (white opaque bottom, tissue culture treated)
- Cell Culture media (DMEM + 10% FBS + 10 mM HEPES + 1x Penicillin-Streptomycin) (Corning Cat# 10017CM, Summerlin Cat# SS-100, Corning Cat# 25-060-CI, Corning Cat# 30-002-CI)
- 0.05% Trypsin-EDTA (Corning Cat# 25-052-CI)
- Firefly Luciferase Assay System (Promega #E1500)

NOTE: Catalog numbers correspond to our recommended products. They can be substituted as desired.

*NOTE: The luciferase RVPs contain luciferase from fireflies (*Photinus pyralis*), therefore a reagent compatible with firefly luciferase is required. Detection reagents for *Renilla luciferase* will not yield a positive result.*

A. Preparation of 293T-hsACE2 cells (see Note 1)

- Lift a sub-confluent flask of 293T-hsACE2 cells using 0.05% Trypsin-EDTA, stop trypsinization by adding Cell Culture media (twice the volume of the Trypsin-EDTA added), pellet the cells by centrifuging in a conical tube for 5 min at 200 g, and resuspend in Cell Culture media.
- Determine cell density using a hemocytometer or equivalent.

B. Thawing, diluting, and plating RVPs

- Remove an aliquot of frozen RVPs, place in a 37°C water bath for 3 minutes until just thawed, and then immediately place on ice (**see Note 2**).
- Gently mix the RVPs tube by inversion.
- Set up RVP dilutions:
 - For 96 well plate, serially dilute the SARS-CoV-2 RVPs to a final volume of 100 µL in Cell Culture media per well. A suggested plate layout is below, with each sample run in duplicate.
 - For 384 well plate, serially dilute the SARS-CoV-2 RVPs to a final volume of 20 µL in Cell Culture media per well.



C. Infection of 293T-hsACE2 cells

a. Infect cells:

- For a 96 well plate, resuspend cells at a density of 2×10^5 cells/mL. Add 100 μ L of the cell suspension (2×10^4 cells) to each well containing the RVPs and mix gently by pipetting 1-2 times (**see Note 3**).
- For a 384 well plate, resuspend cells at a density of 3×10^5 cells/mL. Add 30 μ L of the cell suspension (9×10^3 cells) to each well containing the RVPs and mix gently by pipetting 1-2 times (**see Note 3**).

b. Place plate in a 5% CO₂, 37°C humidified incubator for 48-72 hours (**see Notes 4 and 5**).

D. Quantification of SARS-CoV-2 RVP infection using Promega Luciferase Assay System

a. All subsequent processing of samples should be done using the Promega Luciferase Assay System as per the manufacturer's instructions. Our recommendations are included below.

b. For 96-well assay plates:

- Centrifuge the plates for 5min at 2000rpm to prevent cell loss.
- Equilibrate the 5X lysis buffer to room temperature and then make a 1X solution of lysis buffer using sterile water.
- Carefully aspirate the supernatant with a multichannel pipette.
- Wash cells by carefully adding 100uL of PBS -/- to each well and gently removing as much of the PBS as possible.
- Add 20uL of the lysis buffer to each well of the 96-well assay plate.
- Add 100uL of Luciferase Assay Reagent to each well of the 96-well assay plate.
- Detect Luminescence in a luminometer immediately after reagent addition.

c. For 384-well assay plates:

- Centrifuge the plates for 5min at 2000rpm to prevent cell loss.
- Equilibrate the 5X lysis buffer to room temperature and then make a 1X solution of lysis buffer using sterile water.
- Aspirate supernatant and add 35uL of PBS -/-.
- Aspirate leaving 15uL of liquid behind.
- Add 5uL of the lysis buffer to each well of a 96-well and pulse plate in a centrifuge up to 1000rpm
- Add 20uL of Luciferase Assay Reagent to each well of a 96-well and pulse plate in a centrifuge up to 1000rpm
- Detect Luminescence in a luminometer immediately after reagent addition.

Note: Devices used to read luminescence will vary in relative light unit values based on their individual detectors and software, but the signal to background and Z' will be comparable across devices. Signal to background and Z' can be calculated using mock infected cells as the negative control value.



- A sample 96-well plate layout is below, with each RVP sample in duplicate at 7 concentrations. The final volume in each well is 200 μ L (100 μ L serially diluted RVPs, plus 100 μ L cells).

[illegible]