

Engineering ZIKV E protein to diminish anti-DENV immune response

Identifying the 7 residues that drive DENV-ZIKV cross-reactivity

Ariadna Grinyo-Escuer¹, Sonya Jacobsen¹, Allison Sheetz¹, L. Joe Stafford¹, Kyle Doolan¹, Diana M. Norden¹, Nolan Frey¹, Frederick W. Holtsberg², M. Javad Aman², Edgar Davidson¹, **Benjamin J. Doranz¹**


¹Integral Molecular; ²IBT Bioservices

The high homology between ZIKV and DENV envelope proteins results in immunological cross-reactivity that can exacerbate disease, complicating ZIKV vaccine development

Identifying the 7 ZIKV E protein residues responsible for cross-reactivity enabled engineering of a protective ZIKV immunogen with diminished anti-DENV immune response


Grinyo-Escuer *et al.*, August 2025 *Cell Reports*

171 anti-DENV MAb epitopes




Epitope mapping of MAbs on DENV envelope proteins. Epitope residues shown (in red) summarize 171 MAbs mapped previously by shotgun mutagenesis.

71 MAbs cross-react with ZIKV




DENV MAb binding was tested against ZIKV prM/E expressed in HEK-293T cells (showing DII MAbs only)



Epitopes of 66 cross-reactive MAbs included at least 1 of 7 E protein residues, establishing these 7 residues as the major determinants of DENV-ZIKV cross-reactivity.


"ZIKVm10" E protein engineered to reduce MAb cross-reactivity



(A) The 7 ZIKV E residues responsible for cross-reactivity were mutated to all other amino acids. ELISA was used to identify changes that support SVP production.


(B) Permissive mutations were combined to make "ZIKVm10", which was incorporated into SVPs and tested for binding to ZIKV and DENV MAbs.

ZIKVm10 SVPs used to immunize mice




- CD-1 mice immunized with SVPs displaying WT or ZIKVm10 E protein
- Mouse sera tested for binding with cells expressing DENV or ZIKV prM/E
- Passive immunization of AG129 mice with sera from mice immunized with WT or ZIKVm10 SVPs, then challenged with live ZIKV

ZIKVm10 immunization reduces cross-reactivity




Mean reactivity of sera from 15 naive CD-1 mice, or immunized with WT or ZIKVm10 SVPs, tested by flow cytometry against cells expressing prM/E from A) ZIKV and B) DENV2 (* $p < 0.05$ vs. naive, ** $p < 0.01$, *** $p < 0.001$).

ZIKVm10 immunization maintains protection from live ZIKV




AG129 mice received passive immunizations of serum from naive mice, mice immunized with WT SVPs or ZIKVm10 SVPs, and were then challenged with live ZIKV. Mice were monitored for survival (A), viral load in serum 4 days post-infection (B), body weight, and health score (not shown).

Sera from mice immunized with ZIKVm10 SVPs show reduced ADE of DENV infectivity




Representative ADE-mediated infectivity for sera from individual mice, immunized with WT ZIKV SVPs or ZIKVm10 SVPs, incubated with DENV2 or ZIKV RVPs and K562 cells, mean $\pm$ range from duplicate assays.




For studying immune responses to flaviviruses, handling live virus can be time-consuming with limited assay format or readout options

Reporter Virus Particles are ready-to-use live virus replacements offering optical readouts and flexibility in assay format



Dengue Reporter Virus Particles (RVPs)



Validated Flavivirus RVPs

- Dengue Virus (DENV) Serotypes 1-4
- Zika Virus (ZIKV)
- Yellow Fever Virus (YFV)
- West Nile Virus (WNV)
- 100+ Custom Strains and Mutants


- Capsid, prM/E expressed with modified replicon
- Luciferase/GFP reporter RNA replaces replicon structural genes
- Replication-incompetent flavivirus particles
- Antigenically equivalent to live virus

Applications

- Antibody neutralization
- Serum screening
- Automated workflows
- High-throughput assays


RVPs are highly flexible and have been used for MAb isolation and to study infectivity, neutralization, enhancement (ADE), binding by ELISA, and biosensor

Dengue RVPs




Mattia *et al.*, 2011 *PLoS One*

Stable after thawing, lot-to-lot consistency, good cryopreservation



Work supported by NIAID contracts HHSN272201400058C and HHSN75N93019C00073

RVPs in analyses of response to a Tetravalent Dengue Vaccine



RVP EC₅₀ titers for DENV-1, -3, and -4 after DENV-2 Ab depletion of sera from 71 individuals vaccinated with tetravalent vaccine TAK-003. Sera from 2 studies (red dots or blue dots), collected 30 days after 2nd vaccination, were depleted with mock antigen or DENV-2 antigen, and neutralization titers measured by Dengue RVP assay.

DeMaso *et al.*, 2022 *J Inf Disease*



- Quality manufacturing for clinical trials: ISO 9001 certified processes
- Rapid response to new pathogens
- Ready-to-use convenience, with protocols and support

- Non-replicative reagents available for use in BSL-2 facilities
- 20+ years expertise in virology
- Custom reagents and assays to meet your needs

Benjamin J. Doranz, Ph.D.
215-966-6061
bdoranz@integralmolecular.com
www.integralmolecular.com

Integral Molecular has 20+ years of experience providing virus particles and related services. We support scientists combating diverse viruses: SARS-CoV-2, chikungunya, Ebola, Marburg, Nipah, and more. Visit our website for more details and full catalog.

Download the Full Paper

